



Hypoxyprobe, Inc.  
121 Middlesex Turnpike  
Burlington, MA 01803 USA  
www.hypoxyprobe.com

## Hypoxyprobe™-11 Kit

### Hypoxyprobe™-11 Kit

- Kit contents:** Solid pimonidazole HCl (Hypoxyprobe™-1) Anti-pimonidazole mouse IgG1 monoclonal antibody (MAb1) Lot # 4.3.15 (60 micrograms/ml)
- Data Sheet:** Hypoxyprobe™ RED PE Kit (Cat# HP11-x)
- Description:** Mouse IgG1 monoclonal antibody conjugated to Phycoerythrin (PE) fluorophore (HP-Red-PE-MAb) and x mg of solid pimonidazole HCl (x = 100, 200 or 1000 mg).
- Specificity:** Pimonidazole forms adducts with thiol containing proteins at  $pO_2 \leq 10$  mm Hg. HPRedPE-MAB binds to pimonidazole adducts in hypoxic cells in tissues and in culture. Please see [www.hypoxyprobe.com](http://www.hypoxyprobe.com) for details of the Hypoxyprobe system for detecting cell and tissue hypoxia.
- Format:** Each vial of HP-RedPE-MAB contains anti-pimonidazole antibody dissolved in 200 microliters of stabilized buffer at a concentration of approximately 0.50 mg/mL.
- Product Type:** Fluorophor conjugated IgG1 mouse monoclonal antibody (clone 4.3.11.3)
- Fluorophore- Protein Ratio:**  $\geq 1$  APC moiety per IgG1 molecule.
- Fluorescence:** Excitation 495nm and 565nm; Emission max 573 nm. Phycoerythrin has strong red fluorescence. Product Details
- Applications:** Flow cytometry on isolated cells and immunofluorescence on frozen tissue sections. HP-PE-MAB is an alternative to FITC-labeled anti-pimonidazole mouse monoclonal antibody1. 1/20-1/200 dilution is suggested but users should optimize in their system using appropriate negative and positive controls.
- Target Species:** All.
- Product Form:** Purified IgG1 prepared by affinity chromatography
- Buffer:** Proprietary formulation containing stabilizers. 1.