



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com

Hypoxyprobe™-12 Antibody

Hypoxyprobe™-12 Antibody

- Contents: Anti-pimonidazole Rabbit PAb27 High Titer monoclonal antibody Lot # 10-14-15
- Data Sheet: Hypoxyprobe™ PAb 27 High Titer Antisera (Cat# HP12-x)
- Description: Rabbit PAb 27 monoclonal antibody
- Specificity: Pimonidazole forms adducts with thiol containing proteins at $pO_2 \leq 10$ mm Hg. HP-PAb27 binds to pimonidazole adducts in hypoxic cells in tissues and in culture. Please see www.hypoxyprobe.com for details of the Hypoxyprobe system for detecting cell and tissue hypoxia.
- Format: Each vial of HP-PAb27 contains purified high titer anti-pimonidazole rabbit monoclonal antibody dissolved in 200 microliters of stabilized buffer.
- Product Type: Rabbit monoclonal anti-pimonidazole PAb27 (high titer).
- Target Species: All.
- Product Form: Purified rabbit IgG1 prepared by affinity chromatography
- Buffer: Proprietary formulation containing stabilizers. 1.