



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com

Hypoxyprobe™ 25 Concentrated Mab-1 Antibody

(HPI Catalog # HP 25 Concentrated MAb-1)

- Description: Concentrated Mouse IgG1 monoclonal antibody
- Specificity: Pimonidazole forms adducts with thiol containing tissue proteins at $pO_2 \leq 10\text{mm Hg}$. HP 25 MAb-1 binds to pimonidazole adducts in hypoxic cells in mammalian tissue and tissue culture.
- Format: Each vial of HP-25-MAb-1 contains the antibody dissolved in 250 microliters of PBS solution at a concentration of approximately 1ug/mL.
- Applications: Secondary anti mouse antibody needed. IHC, IF, ELISA, flow cytometry, MAA, Western Blotting.
- Target Species: All mammalian species
- Buffer: Phosphate buffered saline solution