



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com

Data Sheet: Hypoxyprobe™ Anti-FITC HRP
(HPI Catalog # anti-FITC-HRP)

Description:	Rabbit anti-FITC polyclonal conjugated with horseradish peroxidase secondary antibody
Specificity:	Binds to FITC-conjugated primary antibodies
Format:	0.5 mg/mL solution
Applications:	IHC, Western Blotting, flow cytometry, ELISA.
Target Species:	All species
Product Form:	Purified IgG conjugated to HRP
Buffer:	Phosphate buffered saline with 1% BSA, 0.02% Thiomersal
Storage:	+4°C