



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com

Hypoxyprobe™ Anti-FITC HRP Antibody

(HPI Catalog # anti-FITC-HRP)

<u>Description:</u>	Rabbit anti-FITC polyclonal conjugated with horseradish peroxidase secondary antibody
<u>Specificity:</u>	Binds to FITC-conjugated primary antibodies
<u>Format:</u>	0.5 mg/mL solution
<u>Applications:</u>	IHC, Western Blotting, flow cytometry, ELISA.
<u>Target Species:</u>	All species
<u>Product Form:</u>	Purified IgG conjugated to HRP
<u>Buffer:</u>	Phosphate buffered saline with 1% BSA, 0.02% Thiomersal
<u>Storage:</u>	+4°C