



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com



Data Sheet: Hypoxyprobe™ anti Rabbit HRP Secondary Antibody
(HPI Catalog # antirabbithrp)

Description: Goat Anti Rabbit IgG1 polyclonal antibody linked to horse radish peroxidase (HRP).

Specificity: Pimonidazole forms adducts with thiol containing tissue proteins at $pO_2 \leq 10$ mm Hg. HP FITC MAb-1 binds to pimonidazole adducts in hypoxic cells in tissue and tissue culture. Anti PAb HRP binds to HP FITC MAb-1. A chromogen (not supplied) is added to activate the HRP group and stain hypoxic areas of tissue brown.

Format: Each vial of anti -PAb-HRP contains the antibody dissolved in 200 microliters of stabilized diluent at a concentration of approximately 0.5 mg/mL.

Stability: 4 weeks at RT
➤ 2 years at 4 Celcius

Suggested Dilutions: Immunohistology - Paraffin not determined
Immunohistology – Frozen 1/20 to 1/200
Immunohistology – Resin not determined
ELISA 1/200 to 1/1000
Western Blotting 1/20 to 1/200
Flow Cytometry: not determined

Applications:

Target Species: All species.

Product Form: Purified Rabbit IgG1 prepared by affinity chromatography

Buffer: Proprietary formulation containing stabilizers.