



Hypoxyprobe, Inc.
121 Middlesex Turnpike
Burlington, MA 01803 USA
www.hypoxyprobe.com

Hypoxyprobe™ PAb2627(R2710-B) Antibody

(HPI Catalog # R2710-B)

- Description: Rabbit IgG1 polyclonal antipimonidazole antibody
- Specificity: Pimonidazole forms adducts with thiol containing tissue proteins at $pO_2 \leq 10$ mm Hg. HP PAb2627(R2710-B) binds to pimonidazole adducts in hypoxic cells in tissue and tissue culture.
- Format: Each vial of HP-PAb2627(R2710-B) contains the affinity purified antibody dissolved in 200 microliters of stabilized diluent.
- Applications: IHC, immunofluorescence, Western Blotting, flow cytometry, ELISA. Secondary antirabbit antibody needed for all these applications.
- Target Species: All species
- Product Form: Purified Rabbit IgG1 prepared by affinity chromatography
- Buffer: Proprietary formulation containing stabilizers.