



Hypoxyprobe, Inc.  
121 Middlesex Turnpike  
Burlington, MA 01803 USA  
[www.hypoxyprobe.com](http://www.hypoxyprobe.com)

## Hypoxyprobe™ PAb27 Antibody

(HPI Catalog # HP PAb27)

- Description: Rabbit IgG1 polyclonal higher titer anti-pimonidazole antibody.
- Specificity: Pimonidazole forms adducts with thiol containing tissue proteins at  $pO_2 \leq 10$ mm Hg. HP PAb27 binds to pimonidazole adducts in hypoxic cells in tissue and tissue culture.
- Format: Each vial of HP-PAb27 contains the twice affinity purified rabbit polyclonal anti pimonidazole antibody dissolved in 200 microliters of stabilized diluent.
- Applications: Secondary anti rabbit antibody required. IHC, IF, ELISA, MAA, or Western Blotting,
- Target Species: All species
- Product Form: Purified Rabbit polyclonal IgG1 prepared by passing it through two affinity purification chromatography columns.
- Buffer: Proprietary formulation containing stabilizers.